SHAPE AND ROLL:

COSMETIC COVER

Jaipur Foot

- Easy to manufacture
- Gives patients the ability to do many of the things they could with an actual foot
 - Running, climbing trees, etc
- Uses simple materials, and coloring can be made to match the native skin color

Photo of Jaipur Foot removed due to copyright restrictions.

Problems with Jaipur Foot

- Lacks toe support
 - Leads to shortened stride
- Use of different materials inside foot can cause deterioration
- Heavy: current design is roughly 800g

Goals of the Shape & Roll Foot

- Incorporate a foot design to resemble the dorsiflexion of a biological foot/ankle in stride
- Reduce the weight of the foot while maintaining the same durability
- Allow for easier squatting

Current S&R Capabilities

Photo of Shape and Roll Foot from Prosthetics Research Laboratory and Rehabilitation Engineering Research Program at Northwestern University Feinberg School of Medicine has been removed due to copyright restrictions. Closely resembles the roll-over contour of a normal stride Far lighter than the Jaipur Foot Low in cost

Goals for the Shape & Roll Project

- 1. Create an ankle design to allow for extreme dorsiflexion while squatting
 - Edward Sung
- 2. Find a way to make a cosmetic cover for the current S&R Foot
 - Parhys Napier and Nicholas Torgerson

Ankle Joint Design



~ 25 degrees range of motion
Material cost < \$20
Labor cost > \$50

Process:

- Machine grooves and angled sides Sand/grind for aesthetics **Future Plan** Decide on actuation:
 - Wire
 - Push-pull rod
- Work on interface



Goal 2: Improving the Cosmetic Cover

 Create a mold around the S&R core and form the cosmetic cover directly on top of it (S&R-Jaipur Hybrid)

 Find a lightweight, flexible and durable material for the cosmetic cover

Current Problems

Previous attempts at S&R feet with cosmetic covers have still been heavy

 One cosmetic shell design was tested by Cornell, but eventually wore out in a couple of months

Methods of Covering S&R

- Make a cosmetic sleeve, and slide the Shape & Roll into it
 - Fill the empty spaces with a foam-like substance
- Create a hybrid Jaipur Foot using the Shape-and-Roll as the core
 - The cosmetic cover is molded directly onto the S&R using a process similar to that of the Jaipur foot

Cosmetic Shell

Advantages:

- Lightweight
- Rapid process
- Inexpensive

Photo of cosmetic shell for a prosthetic foot removed due to copyright restrictions.

Disadvantages:

- Requires filling due to loose-fit
- Not waterproof
- Depending on the material, sleeve can be abrasive on prosthetic
- Not as rigid as hybrid
- Would not allow for a good bond between the keel and the mold

S&R-Jaipur Hybrid

Advantages:

- Completely watertight
- More durable than a shell
- Less chance for wear due to parts rubbing against one another
- One solid foot rather than multiple pieces

Disadvantages:

- Heavy
- Does not allow for squatting

Photos of hybrid Shape & Roll - Jaipur Foot removed due to copyright restrictions.

Hybrid Foot Approach

Materials in S&R-Jaipur Hybrid

- Rubber (vulcanized, natural, crepe) density: 1100 kg/m³
- Teflon density: 2200 kg/m^3
- Wooden block density: ~ 650 kg/ m^3
- Want to use a similar approach, but use alternate materials for a lighter foot
- Combined with ankle mechanism, the foot will allow for dorsiflexion

Materials – First Prototype

Orethane 60

- Durable, yet also flexible
- Can buy in bulk
- Not very toxic
- ~ \$105/gallon (Smooth-On, Inc.)
- Fiber Reinforcement
 - By adding glass fibers to the Urethane 60, we hope to reinforce the foot and make it more durable
 - ~ \$9/1'x1' .01" sheet (McMaster-Carr)

Additional Materials

- Composite Materials
 - Ex: Carbon fiber
 - Pro: Lightweight, extremely durable
 - Con: Potentially damaging to S&R foot
- Synthetic Rubbers
 - Ex: Neoprene
 - Pros: durability over wide temperature ranges; water resistant
 - Cons: price, local availability
- Microcellular Rubber (MCR)
 - Pros: low density; lightweight; soft; high impact absorption
 - Cons: lacks hardness, susceptible to abrasion

Modification for New Ankle Design

- Use a more flexible material to allow full dorsiflexion capabilities during squatting.
- Looking to use a something similar to a sleeve to cover the ankle.
- Sleeve will be integrated into the mold to make one piece.
- Needs to be waterproof
 - Possible material: urethane 40 (~\$100/ gallon)

Molding Modification

- Current methods use rubber blocks to fill in gaps around the S&R core
- We plan to use a lighter "filler" to reduce weight
 - Foam similar to the kind used in the toes
 - Microcellular Rubber (MCR)
 - Polymers (like polyurethane)
 - Cork
- Sonding between keel and mold
 - Received feedback on the importance of bonding the keel to the mold (decreases movement and wear)
 - Plan to either use strong adhesive or integrate the keel into the mold as one piece.

Plan of Action

- We hope to order materials and begin prototyping the week after spring break (3/29)
 - Begin testing (durability, amount of flexion, etc.) Collect data/analyze (week of 4/5)
 - Modify the design and any failed materials with alternate materials/improved design (week of 4/12)
 - Make new prototype(s) based off of success of previous prototype and resume testing (finish tentatively week of 4/26)
 - Integrate ankle mechanism with cosmetic mold; test with ankle cover (late April, early May)
- Finish a final prototype by the end of the semester.

MIT OpenCourseWare http://ocw.mit.edu

EC.722 Special Topics at Edgerton Center:Developing World Prosthetics Spring 2010

For information about citing these materials or our Terms of Use, visit: http://ocw.mit.edu/terms.